THRIVE ORTHOPEDICS

2024

PRODUCT CATALOG



Supporting clinicians and their patients is our highest priority. This is why we do not sell our AFOs directly to patients online. This is why we have no plans to open clinics to compete against our customers. This is why we stand behind our two-year warranty.

We only succeed if you THRIVE.



Service

We treat every patient like family and every customer like a trusted partner.

Warranty

Our two-year warranty is not just lip service. We do not challenge our clinical partners when they make a claim. Instead, we use these interactions as an opportunity to learn, improve, and build trust.

Design

The purpose of our intentional aesthetic is to inspire patient confidence. We select all materials based on what we would want for ourselves or a family member.

Range

One style of AFO is not right for everyone. We are diligently working towards building the widest range of unique AFO designs on the market.

Values

Thrive is owned and operated by its founders. This has given us the freedom to put our customers first.

We would be honored to have the opportunity to earn your business.

Joe & Mike Co-Founders

Carbon Fiber AFOs

Optimally designed for patients with dexterity issues, Thrive's AFOs were all made with the end user in mind. This results in highly functional products that inspire patient compliance through superior aesthetics and quality.

We know that one AFO is not right for every patient. This is why we are developing the widest range of unique, off-the-shelf carbon fiber AFO designs.







UFO-51 Universal AFO

The UFO-51's patent-pending features provide clinicians with fast assembly and unmatched precision.

Flexible Cuff

heat-moldable calf cuff allows for a wider range of anatomical sizes

Patent Pending Track System

allows for 32 fast, precise height adjustments

Lateral Strut

position uniquely takes advantage of the footplate variability

Trimmable Footplate

footplate adjustments range from a US Women's 5 to Men's 14





PART NUMBER	SIZE	FOOT PLATE LENGTH	BRACE HEIGHT
UFOLT	Universal Left	0.75" 10"	11 5" 10"
UFORT	Universal Right	8.75" - 12"	11.5" - 16"

Carbon Fiber AFOs

F5 AFO

Features a posterior calf cuff with a lateral, posterior, spiral strut design that uniquely avoids contact with the lateral malleolus.

PART NUMBER	SIZE	FOOTPLATE LENGTH	BRACE HEIGHT
SXSLT OR SXSRT	XS	8 7/8"	12"
SSLT OR SSRT	S	9 5/8"	13"
SMLT OR SMRT	М	10 1/2"	14"
SLLT OR SLRT	L	11 1/4"	15"
SXLLT OR SXLRT	XL	12"	16"



F5 Medial AFO

Features a posterior calf cuff with a medial, spiral strut design that is ideal for patients with ankle inversion or supination.

PART NUMBER	SIZE	FOOTPLATE LENGTH	BRACE HEIGHT
MSXLT OR XXSRT	XS	8 7/8"	12"
MSLT OR MSRT	S	9 5/8"	13"
MMLT OR MMRT	М	10 1/2"	14"
MLLT OR MLRT	L	11 1/4"	15"
MXLLT OR MXLRT	XL	12"	16"



Carbon Fiber AFOs

F3 AFO

Features an anterior shell with a lateral strut design that falls distal to the lateral malleolus.

PART NUMBER	SIZE	FOOTPLATE LENGTH	BRACE HEIGHT
AXSLT OR AXSRT	XS	9"	13"
ASLT OR ASRT	S	9 3/4"	14"
AMLT OR AMRT	М	10 3/8"	15"
ALLT OR ALRT	L	11"	16"
AXLLT OR AXLRT	XL	11 3/4"	17"



F3 Hero Pediatric AFO

Flexibility and shape specifically designed for the pediatric population, this AFO features an anterior shell with a lateral strut design that falls distal to the lateral malleolus.

PART NUMBER	SIZE	FOOTPLATE LENGTH	BRACE HEIGHT
PASLT OR PASRT	S	6 1/4"	8 5/8"
PAMLT OR PAMRT	М	7"	10 1/8"
PALLT OR PALRT	L	7 7/8"	11 5/8"
PAXLLT OR PAXLRT	XL	8 1/4"	12 3/8"



Equinus Brace

One hour per day treatment for plantar fasciitis and equinus deformity. The Equinus Brace™ is the only orthosis on the market that telescopes above the knee to stretch both the gastrocnemius and soleus muscles. Ankle joints allow for a progressive, measurable treatment protocol. The toe wedge engages the windlass mechanism, providing an optimal stretch through the ankle joint while also stretching the plantar fascia.



"We're using them a lot for our Achilles tendinitis patients we'll often have them in a soft cast and cam walker until pain-free, then we'll transition them into an ankle brace supplemented by daily active stretching and use of the Equinus Brace. Our patients are giving us positive feedback and we're seeing noticeable increases in their flexibility and range of motion after using the brace, so we are very pleased."

James S., DPM

Equinus Brace

	NIGHT SPLINT	EQUINUS BRACE	MANUAL STRETCHING
Locks Knee into Extension	×	~	✓
Controllable Stretch	×	~	×
Not Technique Dependent	×	~	×
1 Hour Per day	×	~	✓
Ease of Use	×	~	X

Patented "Above the Knee Extension"

the only brace that extends above the knee maintaining full extension while dorsiflexing the foot.

Adjustable Hinges

allow for precise, gradual, and safe stretching while providing consistent and measurable results.

Patented "Toe Wedge"

Provides an additional stretch directly to the plantar fascia and correctly positions the foot and ankle for an optimal stretch.



PART NUMBER	SIZE	DESCRIPTION
IQ1004	US TODDLER 8 - US LITTLE KIDS 2	UNIVERSAL PEDIATRIC
IQ1001	US W6 - US M11	UNIVERSAL ADULT
IQ1003	US M11 - M 16	UNIVERSAL ADULT XL



Our foot plates are 100% carbon fiber, providing superior strength through a continuous weave with true dynamic feedback. The flat foot plates are conveniently sold in 12 shoe sizes, saving clinicians significant adjustment time.





*Insoles sold individually

	FULL LENGTH	FULL LENGTH	MORTON'S EXT	MORTON'S EXT			
SHOE SIZE	PART NO. 1.3 MM	PART NO. 1.0 MM	PART NO. 1.3 MM	PART NO. 1.0 MM	LENGTH	HEEL WIDTH	MET HEAD WIDTH
W5	FPW5	FPFW5	RMEW5	FMEW5	8.4" 213mm	1.8" 45mm	1.9" 73mm
W6	FPW6	FPFW6	RMEW6	FMEW6	8.8" 223mm	1.9" 47mm	3" 75m
W7 M6	FPM6	FPFM6	RMEM6	FMEM6	9.2" 233mm	1.9" 48mm	3.1" 78mm
W8 M7	FPM7	FPFM7	RMEM7	FMEM7	9.6" 243mm	2" 50mm	3.2" 81mm
W9 M8	FPM8	FPFM8	RMEM8	FMEM8	10" 253mm	2" 52mm	3.3" 84mm
W10 M9	FPM9	FPFM9	RMEM9	FMEM9	10.4" 263mm	2.1 54mm	3.4" 87mm
W11 M10	FPM10	FPFM10	RMEM10	FMEM10	10.7" 273mm	2.2" 55mm	3.5" 89mm
W12 M11	FPM11	FPFM11	RMEM11	FMEM11	11.1" 283mm	2.2" 56mm	3.5" 90mm
W13 M12	FPM12	FPFM12	RMEM12	FMEM12	11.5" 293mm	2.2" 57mm	3.6" 91mm
M13	FPM13	FPFM13	RMEM13	FMEM13	11.9" 302mm	2.2" 57mm	3.6" 92mm
M14	FPM14	FPFM14	RMEM14	FMEM14	12.3" 312mm	2.3" 58mm	3.7" 94mm
M15	FPM15	FPFM15	RMEM15	FMEM15	12.7" 322mm	2.3" 58mm	3.7" 94mm





*Insoles sold individually

SIZE	SHOE SIZE RANGE	PART NO. 1.0 MM	PART NO. 1.3 MM	LENGTH	HEEL WIDTH	MET HEAD WIDTH
S	Up to US W6	FCCIS	CCIS	9" 228mm	2" 51mm	2.4" 71mm
M	US W6.5 to US M7	FCCIM	CCIM	10"253mm	2.1" 53mm	3.1" 78 mm
L	US M7.5 to US M10	FCCIL	CCIL	11" 280mm	2.3" 59mm	3.4" 86mm
XL	US M10.5 - US M13.5	FCCIXL	CCIXL	12.1" 308mm	2.7" 69mm	3.6" 91mm

Carbon FLX

*Insoles sold as a pair

SHOE SIZE	PART NUMBER	LENGTH	HEEL WIDTH	MET HEAD WIDTH
W7 M6	FLXW7M6	9.3" 237mm	2.4" 61.3mm	3.4" 86mm
W8 M7	FLXW8M7	9.7" 247mm	2.5" 63.6mm	3.5" 89mm
W9 M8	FLXW9M8	10.1" 256.5mm	2.6" 65.9mm	3.6" 91mm
W10 M9	FLXW10M9	10.5" 267mm	2.6" 67mm	3.7" 95mm
W11 M10	FLXW11M10	10.9" 277mm	2.7" 69mm	3.9" 98mm
W12 M11	FLXW12M11	11.3" 286.5mm	2.8" 71.5mm	4" 101mm
W13 M12	FLXW13M12	11.8" 299mm	2.9" 72.4mm	4" 102.5mm
M13	FLXW14M13	12.2" 309.3mm	2.9" 74mm	4.1" 103.9mm
M14	FLXM14	12.5" 317.5mm	3" 74.8mm	4.3" 108.5mm
M15	FLXM15	12.9" 327.5mm	3.1" 78mm	4.3" 109.5mm



Ankle Braces

The Jetlace features a unique rear-entry design with removable medial and lateral ankle supports, ensuring easy donning for patients with dexterity issues. Both models utilize a ballistic nylon shell coupled with our moisture-wicking internal padding. Every Thrive ankle brace includes a complimentary wash bag for easy cleaning.

Ankle Braces

Mach8[™]

PART NUMBER	SIZE	SHOE SIZE
MACH-XS	XS	M 6-7 W 7-8
MACH-S	S	M 7.5-9 W 8.5-10
MACH-M	М	M 9.5-11 W 10.5-12
MACH-L	L	M 11.5-13 W 12.5-14
MACH-XL	XL	M 13.5 -15 W 14.5+



Jetlace[™]

PART NUMBER	SIZE	SHOE SIZE
JET-XS	XS	M 3-5 W 4-6
JET-S	S	M 6-7 W 6.5-9
JET-M	М	M 7.5-10 W 9.5-11
JET-L	L	M 10.5-12 W 11.5-13
JET-XL	XL	M 12.5 -14 W 13.5-15



Accessories

Magnetic Strapping System

Designed for low maintenance longevity, this upgrade is a great option for reducing the need to replace velcro. In addition, it is a great option for patients with dexterity issues.



Symmetry T-Strap[™]

Add to any Thrive AFO to stabilize mild to moderate ankle instability.



Accessories

CoolFit AFO Socks™

Padded Coolmax $^{\text{TM}}$ provides cooling comfort and protection for AFO users.




X-Fit Calf Sleeve

Provides compression of 20-30 mmHG and provides comfort and protection for AFO users, while allowing them to use their preferred choice of socks.



Testimonials



"Resilience is accepting your new reality"

Despite drop foot, I was going to run again and complete this goal. 7 months postpartum, I am back running in races weekly. It's different now but so much more rewarding. I run with foot drop, while pushing my 7-month-old baby girl in her stroller.

AFO User - Shelby G.

"They are 100% behind their product"

"Any minute there is an issue they are 100% behind their product. Very nice to see these days when nobody gives a s***."

Customer Service - Erik S., CPO

"The fit and finish is superb"

"The fit and finish of the Thrive AFO is superb. With the use of the magnetic closures the patients with limited hand dexterity are assisted both in donning and doffing. A great product!

Magnetic Straps - Michael J., CPO

"We tried injections, physical therapy, night braces etc etc. Nothing worked. Your brace is helping him."

"I just want to tell you how excited I am about your brace. One of my patients, whom I have been treating for plantar fasciitis for months, is now for the first time showing signs of improvement. We tried injections, physical therapy, night braces etc etc. Nothing worked. Your brace is helping him."

Equinus Brace - Elliot U., DPM

Testimonials

"As a professional in design..., I appreciate the sleek design."

As a professional in the design and construction industry, I wanted to get back to traveling for work and performing site visits. I am also a fan of walking outdoors, especially on the greenways in my city. I am extremely happy with my F5 Spiral AFO. The padding is comfortable and I appreciate the sleek design.

AFO User - Lansing D.

"It is a pleasure to work with such a competent company"

"Working with Thrive Orthopedics with regards to placing the order and receiving it within a short period is imperative to the welfare of our patients. It is a pleasure to work with such a competent company and the AFO is just as great!"

Turnaround Time - Charlene G., Office Manager

"Your brace gives me hope"

"I'm 36 years old, have had 6 ankle surgeries, have little cartilage left, and struggle to even take long walks without limping around for a week. Riding my bike and using the erg machine has even become a challenge. Your brace gives me hope that I can run around with my son who isn't even a year old and helps me push back a full ankle replacement for as long as possible. I commend you all for everything you do. Keep kicking ass, thank you!!"

AFO User - John N.



MORE SUPPORT

Decision Guide

Does the patient need additional support for any of the following conditions?

- Varus or Valgus Knee Instability
- Quadriceps Weakness
- Foot, Ankle AND Knee Instability
- Triplanar Knee and Ankle Instability

Does the patient have knee flexion

instability or a partial foot amputation?

Does the patient have ankle eversion/

pronation or a medial bony prominence?

CONSIDER...

Hercules KAFO

Suggested L2036 Custom fabricated



CONSIDER...

Magnetic Available



F3 Anterior AFO

PDAC Approved L1932



CONSIDER...

F5 Spiral AFO

PDAC Approved L1951 Magnetic Available



Does the patient have ankle inversion/supination or lateral bony prominence?

CONSIDER...

F5 Medial AFO

PDAC Approved L1951 Magnetic Available



KNEE HYPEREXTENSION


CONSIDER...

UFO-51 Universal AFO

PDAC Pending L1951 Magnetic Available



Will the patient need incremental height adjustments for the treatment plan?





Scan the QR Code, Email or Call to order.



sales@thriveorthopedics.com

TEL: 484.442.0494

FAX: 267.573.9765

